

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human RPUSD1 protein



## Basic Information

**Catalog Number:**

Ag22314

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET30a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MEPGSVENLSIVYRSRDFLVVNKHWVDVRI~~DSKAWRET~~  
LTLQKQLRYRFP~~ELADPDTCYGF~~R~~FCHQLDFSTSGAL~~  
CVALNKAAAGSAYRC~~FKERRVT~~KAYLALLRGHIQESR  
VTISHAIGRNSTEGRAHTMCIEG~~SQGCENPKPSLTDL~~  
VVLEHG~~LYAGDPVSKVLLKPLTGRTHQLRVHCSALGH~~  
PVVGD~~LT~~YGEVSGREDRPF~~RMMLHAFYLRIPTDTECV~~  
EVCTPDP~~FLPSLDACWSPHTLLQSLDQLVQALRATPD~~  
PDPE~~DRGPRPGSPSALLPGRPPPPPTKPPETEAQ~~RG  
PCLQWLSEW~~TLEPDS~~

(1-312 aa encoded by BC011783)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W:  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.