

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ULK2 protein



## Basic Information

Catalog Number:

Ag3695

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MEVVGDFEYSKRDLVGHGAFAVVFRGRHRQKTDWE  
VAIKSINKKNLSKSQILLGKEIKILKELQHENIVALYDVQ  
ELPNSVFLVMEYCNGGDLADYLQAKGTLS EDTIRVFL  
HQIAAAMRILHSGIHRDLKPQNILLSYANRRKSSVSG  
IRIKIADFGFARYLHSNMMAATLCGSPMYMAPEVIMS  
QHYDAKADLWSIGTVIYQCCLVGKPPFQANS PQDLRM  
FYEKNRSLMPSIPRETSPYLANLLGLLQRNQKDRMDF  
EAFFSHPFLEQGPV

(1-276 aa encoded by BC034988)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.