

For Research Use Only

Biotin Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody

Catalog Number: Biotin-65607



Basic Information

Catalog Number:

Biotin-65607

Size:

100ug, 500 µg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

W3/13

Applications

Tested Applications:

FC

Species Specificity:

rat

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 rat splenocytes were surface stained with 0.25 μ g Biotin Anti-Rat CD43 (W3/13) Mouse IgG2a RecAb (Biotin-65607, Clone: W3/13) and CoraLite®488-conjugated streptavidin, or unstained. Cells were not fixed.



1×10^6 rat splenocytes were surface stained with 0.25 μ g Biotin Anti-Rat CD43 (W3/13) Mouse IgG2a RecAb (Biotin-65607, Clone: W3/13) and CoraLite®488-conjugated streptavidin (red), or unstained (blue). Cells were not fixed.