

For Research Use Only

CoraLite® Plus 405 Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody



Catalog Number: **CL405-65607**

Basic Information

Catalog Number:

CL405-65607

Size:

100ug , 500 µg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

W3/13

Excitation/Emission maxima
wavelengths:

399 nm / 422 nm

Applications

Tested Applications:

FC

Species Specificity:

rat

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 rat splenocytes were surface stained with 0.25 ug CoraLite® Plus 405 Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody (CL405-65607, Clone:W3/13), or CoraLite® Plus 405 Mouse IgG2a Isotype Control (C1.18.4) (CL405-65208, Clone: C1.18.4). Cells were not fixed.



1×10^6 rat splenocytes were surface stained with 0.25 ug CoraLite® Plus 405 Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody (CL405-65607, Clone:W3/13) (red), or CoraLite® Plus 405 Mouse IgG2a Isotype Control (C1.18.4) (CL405-65208, Clone: C1.18.4) (blue). Cells were not fixed.