

For Research Use Only

CoraLite® Plus 488-conjugated RPE Polyclonal antibody



Catalog Number: **CL488-12168**

Basic Information

Catalog Number:

CL488-12168

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2810

GenBank Accession Number:

BC016764

GeneID (NCBI):

6120

UNIPROT ID:

Q96AT9

Full Name:

ribulose-5-phosphate-3-epimerase

Calculated MW:

228 aa, 25 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

