

For Research Use Only

CoraLite® Plus 488-conjugated TRBP Polyclonal antibody

Catalog Number: CL488-15753

Featured Product



Basic Information

Catalog Number:

CL488-15753

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8255

GenBank Accession Number:

BC005860

GeneID (NCBI):

6895

UNIPROT ID:

Q15633

Full Name:

TAR (HIV-1) RNA binding protein 2

Calculated MW:

39 kDa

Observed MW:

43-46 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:50-1:500

**Excitation/Emission maxima
wavelengths:**

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC, FC (Intra)

Species Specificity:

human, mouse, rat

Positive Controls:

IF/ICC : HeLa cells,

Background Information

TARBP2 is a component of the RISC loading complex (RLC), also known as the micro-RNA (miRNA) loading complex (miRLC). Within the complex, TARBP2 and DICER1 are required to process precursor miRNAs (pre-miRNA) to mature miRNA and then load them onto EIF2C2/AGO2. It may have a role in the production of short interfering RNAs (siRNAs) from double-stranded RNA (dsRNA) by DICER1. Otherwise, it binds to the HIV-1 TAR RNA that is located in the long terminal repeat (LTR) of HIV-1, and stimulates translation of TAR-containing RNAs. The endogenous TRBP that was present in HEK293 cells, with apparent molecular masses of 43-46 kDa, are expressed in human cells. (PMID: 16142218)

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using Coralite® Plus 488 TRBP antibody (CL488-15753) at dilution of 1:200, CL594-Phalloidin (red).



1x10⁶ HeLa cells were intracellularly stained with 0.8 ug Coralite® Plus 488 TRBP Polyclonal Antibody (CL488-15753)(red), or 0.8 ug Coralite® Plus 488-conjugated Rabbit IgG control Rabbit PolyAb (CL488-30000) (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).