

For Research Use Only

CoraLite® Plus 488-conjugated TRIM15 Recombinant antibody

Catalog Number:CL488-83582-5



Basic Information

Catalog Number: CL488-83582-5	GenBank Accession Number: BC038585	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 89870	CloneNo.: 240565C7
Source: Rabbit	UNIPROT ID: Q9C019	Recommended Dilutions: IF/ICC 1:50-1:500
Isotype: IgG	Full Name: tripartite motif-containing 15	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG4533	Calculated MW: 465 aa, 52 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : HeLa cells,
Species Specificity: human	

Background Information

TRIM15 (Tripartite motif-containing protein 15) is also named as RNF93, ZNF178, ZNFB7 and belongs to the TRIM/RBCC family. It acts as an ubiquitin-protein ligase and contributes to the endogenous restriction of retroviruses in cells.

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CoraLite® Plus 488 TRIM15 antibody (CL488-83582-5, Clone: 240565C7) at dilution of 1:200, CL594-Phalloidin (red).