

For Research Use Only

# CoraLite® Plus 488-conjugated MARS Recombinant antibody

Catalog Number:CL488-83690-5



## Basic Information

<b>Catalog Number:</b> CL488-83690-5	<b>GenBank Accession Number:</b> BC002384	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 4141	<b>CloneNo.:</b> 240648G3
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P56192	<b>Recommended Dilutions:</b> IF/ICC 1:50-1:500
<b>Isotype:</b> IgG	<b>Full Name:</b> methionyl-tRNA synthetase	<b>Excitation/Emission maxima wavelengths:</b> 493 nm / 522 nm
<b>Immunogen Catalog Number:</b> AG6619	<b>Calculated MW:</b> 101 kDa	
	<b>Observed MW:</b> 101 kDa	

## Applications

<b>Tested Applications:</b> IF/ICC	<b>Positive Controls:</b> IF/ICC : HepG2 cells,
<b>Species Specificity:</b> human, mouse, rat	

## Background Information

The methionyl-tRNA synthetase (MARS) gene encodes cytoplasmic methionyl-tRNA synthetase (MetRS) responsible for catalyzing the ligation of methionine to tRNA. MetRS belongs to a family of aminoacyl-tRNA synthetases that play critical roles in protein biosynthesis by charging tRNAs with their cognate amino acids. MetRS overexpression was shown to be evident in human colon cancer patients. MetRS may thus be involved in oncogenic transformation. MetRS has been proven as a potential prognostic marker candidate for the clinical prognostic prediction of non-small-cell lung cancer (NSCLC) in patients. MetRS can also be used to detect intracellular oxidative stress and control protein synthesis under oxidative stress (PMID: 30271085, PMID: 34679529, PMID: 32404475). The calculated molecular weight of MetRS is 101 kDa.


## Storage

**Storage:**  
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using CoraLite® Plus 488 MARS antibody (CL488-83690-5, Clone: 240648G3) at dilution of 1:200.