

For Research Use Only

CoraLite® Plus 488-conjugated ANKRD6 Recombinant antibody

Catalog Number: CL488-84210-4



Basic Information

Catalog Number:

CL488-84210-4

GenBank Accession Number:

BC042173

Purification Method:

Protein A purification

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

GeneID (NCBI):

22881

CloneNo.:

241494D11

Source:

Rabbit

UNIPROT ID:

Q9Y2G4

Recommended Dilutions:

IF/ICC 1:50-1:500

Isotype:

IgG

Full Name:

ankyrin repeat domain 6

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Immunogen Catalog Number:

AG19328

Calculated MW:

727 aa, 80 kDa

Applications

Tested Applications:

IF/ICC

Positive Controls:

IF/ICC : A549 cells,

Species Specificity:

human

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed A549 cells using Coralite® Plus 488 ANKRD6 antibody (CL488-84210-4, Clone: 241494D11) at dilution of 1:200.