

For Research Use Only

# CoraLite® Plus 488-conjugated THAP11 Recombinant antibody

Catalog Number:CL488-84373-6



## Basic Information

Catalog Number:

CL488-84373-6

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19342

GenBank Accession Number:

BC012182

GeneID (NCBI):

57215

UNIPROT ID:

Q96EK4

Full Name:

THAP domain containing 11

Calculated MW:

314 aa, 34 kDa

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

241590E5

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse

Positive Controls:

IF/ICC : HeLa cells,

## Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using Coralite® Plus 488 THAP1 antibody (CL488-84373-6, Clone: 241590E5) at dilution of 1:200.