

For Research Use Only

# CoraLite®594 Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody

Catalog Number: CL594-65571



## Basic Information

Catalog Number:

CL594-65571

Size:

100ug , 500 µg/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

MOPC-173

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

## Applications

Tested Applications:

FC

Species Specificity:

n/a

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human PBMCs were unstained (blue) or surface stained with 5 ul CoraLite®594 Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody (CL594-65571, Clone: MOPC-173) (red) or CoraLite® Plus 594 Anti-Human CD3 (Hit3a) (CL594-65112, Clone: Hit3a) (green). Cells were not fixed.