

For Research Use Only

CoraLite®594 Anti-Human CD16 (3G8) Mouse IgG2a Recombinant Antibody

Catalog Number: CL594-65612



Basic Information

Catalog Number:
CL594-65612

Size:
100 tests, 5 ul/test

Source:
Mouse

Isotype:
IgG2a

GenBank Accession Number:
BC017865

GeneID (NCBI):
2214

ENSEMBL Gene ID:
ENSG00000203747

Full Name:
Fc fragment of IgG, low affinity IIIa,
receptor (CD16a)

Calculated MW:
254 aa, 29 kDa

Purification Method:
Protein A purification

CloneNo.:
3G8

Excitation/Emission maxima
wavelengths:
588 nm / 604 nm

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

CD16 is a 50-70-kDa low affinity Fc receptor found on the surface of natural killer cells, neutrophil polymorphonuclear leukocytes, monocytes and macrophages. CD16 mediates antibody-dependent cellular cytotoxicity (ADCC) and other antibody-dependent responses, such as phagocytosis. CD16 has been identified as Fc receptors FcγRIIIa (CD16a) and FcγRIIIb (CD16b), encoded by two nearly identical genes, FCGR3A and the FCGR3B. Clone 3G8 recognizes both the CD16a and CD16b (PMID: 7592758).

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human peripheral blood leukocytes were surface stained with 5 ul CoraLite®594 Anti-Human CD16 (3G8) Mouse IgG2a RecAb (CL594-65612, Clone:3G8), or CoraLite®594 Mouse IgG2a Isotype Control (C1.18.4) (CL594-65208, Clone: C1.18.4). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.



1x10⁶ human peripheral blood leukocytes were surface stained with 5 ul CoraLite®594 Anti-Human CD16 (3G8) Mouse IgG2a RecAb (CL594-65612, Clone:3G8)(red), or CoraLite®594 Mouse IgG2a Isotype Control (C1.18.4) (CL594-65208, Clone: C1.18.4) (blue). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.