

For Research Use Only



CoraLite®594-conjugated Integrin Alpha V Recombinant monoclonal antibody

Catalog Number:CL594-84883-5

Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
CL594-84883-5	BC136442	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
100ul , Concentration: 1000 ug/ml by Nanodrop;	3685	242376D3
Source:	UNIPROT ID:	Recommended Dilutions:
Rabbit	P06756	IF/ICC: 1:50-1:500
Isotype:	Full Name:	Excitation/Emission maxima
IgG	integrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51)	wavelengths: 588 nm / 604 nm
Immunogen Catalog Number:	Calculated MW:	
AG25825	1048 aa, 116 kDa	
	Observed MW:	
	130 kDa	

Applications

Tested Applications:	Positive Controls:
IF/ICC	IF/ICC : A549 cells,
Species Specificity:	
human	

Background Information

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated alpha and beta subunits. Integrin alpha V, also known as CD51, is a type I transmembrane glycoprotein composed of an extracellular domain, a transmembrane region and cytoplasmic domain (PMID: 2430295; 1690718). It can form heterodimers with one of the five beta integrin subunits (beta 1, 3, 5, 6, 8) that recognize ligands containing the sequence Arg-Gly-Asp, such as vitronectin, fibronectin, and osteopontin (PMID: 37087799). Integrin alpha V mainly plays a role in extracellular matrix-mediated cell adhesion and migration, differentiation, wound healing, inflammation, tumorigenesis, proliferation, angiogenesis, and metastasis.


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed A549 cells using CoraLite®594 Integrin Alpha V antibody (CL594-84883-5, Clone: 242376D3) at dilution of 1:200, CL488-Phalloidin (green).