

For Research Use Only

# CoraLite® Plus 647 Anti-Mouse CD48 (HM48-1)



Catalog Number: **CL647-65227**

## Basic Information

Catalog Number:

CL647-65227

Size:

100ug, 0.5 mg/ml

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

BC060977

GeneID (NCBI):

12506

UNIPROT ID:

P18181

Full Name:

CD48 antigen

Purification Method:

Protein A purification

CloneNo.:

HM48-1

Excitation/Emission maxima  
wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface stained with 1 ug CoraLite® Plus 647 Anti-Mouse CD48 (CL647-65227, Clone: HM48-1) (red), or 1 ug Control Antibody. Cells were not fixed.