

For Research Use Only

CoraLite® Plus 647-conjugated Di-Methyl-Histone H3 (Arg2) Recombinant monoclonal antibody

Catalog Number: **CL647-84702**



Basic Information

Catalog Number:

CL647-84702

GenBank Accession Number:

BC066245

Purification Method:

Protein A purification

Size:

100ul , Concentration: 1000 ug/ml by

Nanodrop;

GeneID (NCBI):

8350

CloneNo.:

241896D7

Source:

Rabbit

UNIPROT ID:

P68431

Recommended Dilutions:

IF/ICC: 1:50-1:500

Isotype:

IgG

Full Name:

histone cluster 1, H3a

Excitation/Emission maxima

wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

IF/ICC

Positive Controls:

IF/ICC : HEK-293 cells,

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HEK-293 cells using CoraLite® Plus 647 Di-Methyl-Histone H3 (Arg2) antibody (CL647-84702, Clone: 241896D7) at dilution of 1:200.