

For Research Use Only

CoraLite® Plus 647-conjugated Acetyl-Histone H3 (Lys4) Recombinant antibody

Catalog Number: **CL647-84831**



Basic Information

Catalog Number: CL647-84831	GenBank Accession Number: BC066245	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 8350	CloneNo.: 241882E9
Source: Rabbit	UNIPROT ID: P68431	Recommended Dilutions: IF/ICC 1:50-1:500
Isotype: IgG	Full Name: histone cluster 1, H3a	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
	Observed MW: 15 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : HeLa cells,
Species Specificity: human, mouse	

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 647 Acetyl-Histone H3 (Lys4) antibody (CL647-84831, Clone: 241882E9) at dilution of 1:200.