

For Research Use Only

# CoraLite® Plus 750 Anti-Human CD49b (AK7)



Catalog Number: CL750-65252

## Basic Information

<b>Catalog Number:</b> CL750-65252	<b>GenBank Accession Number:</b> NM_002203	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 100tests , 5 µl/test	<b>GeneID (NCBI):</b> 3673	<b>CloneNo.:</b> AK7
<b>Source:</b> Mouse	<b>Full Name:</b> integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2 receptor)	<b>Excitation/Emission maxima wavelengths:</b> 755 nm / 780 nm
<b>IsoType:</b> IgG1, kappa	<b>Calculated MW:</b> 126 kDa	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Human

## Background Information

CD49b, also known as alpha2 integrin, is a type I transmembrane glycoprotein which forms an adhesion molecule VLA-2 with beta1 integrin. CD49b is a receptor for many (extracellular) matrix and non-matrix molecules (PMID: 23847623). CD49b is expressed on numerous cell types including platelets, activated T cells, B cells, and monocytes. CD49b activation skews NK cell effector function to a phenotype that exhibits elevated IFN-γ production with a concomitant inhibition of motility and cytolytic function (PMID: 32971468, 32762681). The Anti-Human CD49b antibody, clone AK7, is derived from Human platelets.

## Storage


**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD49b (CL750-65252, Clone: AK7) or unstained. Cells were not fixed. Monocytes were gated.