

For Research Use Only

# CoraLite® Plus 750-conjugated MAPKAPK2 Recombinant antibody

Catalog Number:CL750-82929-2



## Basic Information

Catalog Number:

CL750-82929-2

Size:

100ul , Concentration: 1000 µg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34196

GenBank Accession Number:

BC036060

GeneID (NCBI):

9261

UNIPROT ID:

P49137

Full Name:

mitogen-activated protein kinase-  
activated protein kinase 2

Calculated MW:

400 aa, 46 kDa

Observed MW:

42-46 kDa

Purification Method:

Protein A purification

CloneNo.:

230197D6

Recommended Dilutions:

WB 1:500-1:1000

Excitation/Emission maxima  
wavelengths:

755 nm / 780 nm

## Applications

Tested Applications:

WB

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, HeLa cells

## Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with CL750-82929-2 (MAPKAPK2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.