

For Research Use Only

CoraLite® Plus 750 Anti-Human CD324 (E-cadherin) Rabbit Recombinant Antibody



Catalog Number: **CL750-98030**

Basic Information

Catalog Number: CL750-98030	GenBank Accession Number: BC141838	Purification Method: Protein A purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 999	CloneNo.: 240159G8
Source: Rabbit	UNIPROT ID: P12830	Excitation/Emission maxima wavelengths: 755 nm / 780 nm
Isotype: IgG	Full Name: cadherin 1, type 1, E-cadherin (epithelial)	
Immunogen Catalog Number: AG14973	Calculated MW: 882 aa, 97 kDa	

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

Cadherins are a family of transmembrane glycoproteins that mediate calcium-dependent cell-cell adhesion and play an important role in the maintenance of normal tissue architecture. E-cadherin (epithelial cadherin), also known as CD324, CDH1 (cadherin 1) or CAM 120/80, is a classical member of the cadherin superfamily which also include N-, P-, R-, and B-cadherins. E-cadherin is expressed on the cell surface in most epithelial tissues. The extracellular region of E-cadherin establishes calcium-dependent homophilic trans binding, providing specific interaction with adjacent cells, while the cytoplasmic domain is connected to the actin cytoskeleton through the interaction with p120-, α -, β -, and γ -catenin (plakoglobin). E-cadherin is important in the maintenance of the epithelial integrity, and is involved in mechanisms regulating proliferation, differentiation, and survival of epithelial cell. E-cadherin may also play a role in tumorigenesis. It is considered to be an invasion suppressor protein and its loss is an indicator of high tumor aggressiveness.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**

Selected Validation Data



1x10⁶ HT-29 cells (left) or Jurkat cells (right) were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD324 (E-cadherin) Rabbit Recombinant Antibody (CL750-98030, Clone: 240159G8) (red), or Isotype Control (blue). Cells were not fixed.