

For Research Use Only

# Recombinant Human FOLH1 protein (His Tag)



Catalog Number: Eg0216

## Basic Information

**Species:**  
Human

**Purity:**  
>90 %, SDS-PAGE

**Tag:**  
His Tag

## Technical Specifications

**Purity:**

>90 %, SDS-PAGE

**Endotoxin Level:**

<1.0 EU/μg protein, LAL method

**Source:**

HEK293-derived Human FOLH1 protein Lys44-Ala750 (Accession# Q04609) with a His tag at the N-terminus.

**GeneID:**

2346

**Accession:**

Q04609

**Predicted Molecular Mass:**

80.6 kDa

**SDS-PAGE:**

**Formulation:**

Lyophilized from sterile PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

## Biological Activity

Not tested

## Storage and Shipping

**Storage:**

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

**Shipping:**

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

## Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

## Background

## References

## Synonyms

PSMA, PSMA/GCPII, Cell growth-inhibiting gene 27 protein, EC:3.4.17.21, FGCP

### For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

