

For Research Use Only

FITC Plus Anti-Human STRO-1 (STRO-1)



Catalog Number: FITC-65184

1 Publications

Basic Information

Catalog Number:

FITC-65184

Size:

100tests , 5 µl/test

Source:

Mouse

Isotype:

IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Purified by precipitation and chromatography

CloneNo.:

STRO-1

Excitation/Emission maxima wavelengths:

495 nm / 524 nm

Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

Human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Tao Jiang	37441553	J Tissue Eng	FC

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human whole blood cells were surface stained with 5 ul FITC Plus Anti-Human STRO-1 (FITC-65184, Clone: STRO-1) or unstained. Cells were not fixed. CD45dim, CD34+ cells were gated.