

For Research Use Only



MultiPro™ 5CFLX Mouse IgG1 Isotype Control (1F8D3)

Catalog Number: G66360-1-5C

Basic Information

Catalog Number:

G66360-1-5C

Size:

10ug, Concentration: 500ug/mL by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

Full Name:

MultiPro™ 5CFLX Mouse IgG1 Isotype Control (1F8D3)

CloneNo.:

1F8D3

Conjugate:

5CFLX

Full Oligo Sequence:

CGGAGATGTGTATAAGACAGTCTG
AATATGTCGTGCCCATATAAGAAA

Barcode Sequence:

TCTGAATATGTCGTG

Applications

Tested Applications:

Single Cell (Intra)

Species Specificity:

NA

Background Information

This antibody is used as the isotype control of mouse IgG1 (Kappa). It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

Storage

Storage:

2-8°C

Storage Buffer:

PBS with 1mM EDTA and 0.09% sodium azide


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



G66360-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G66360-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.