

For Research Use Only

MultiPro® Enhanced Blocking Reagent

Catalog Number: G900005



Tested Applications

Positive SC detected in [10x Genomics Chromium Gene Expression Flex with Feature Barcodes](#).

Enhanced Blocking Reagent is designed to decrease antibody background staining when performing.

Package


65uL

Recommended dilution

APPLICATION	DILUTION
SC intracellular staining	15uL per 2M cells

Sample-dependent, Check data in validation data gallery.

Validation Data



SC experiment of Human PBMCs using G900005:
Human PBMCs were stimulated with PMA for 4 hours. Cells were blocked with either Enhanced Blocking Reagent or BSA.

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.**