

## HumanKine<sup>®</sup> recombinant human IL-9 protein- GMP grade



Animal Component-Free

Human cell expressed

Tag-Free

Endotoxin Free

### Product Description


Animal-free Recombinant Human IL-9 is expressed from human 293 cells as a glycoprotein monomer with apparent molecular mass of 38 to 48 kDa. IL-9 is produced in a human cell expression system with serum-free, chemically defined media. It stimulates cell proliferation and also prevents apoptosis. IL-9 acts as a regulator for a variety of hematopoietic cells. It enhances the expansion and recruitment of mast cells and eosinophils. This cytokine is greater than 95% pure.

Alternative Names	Cytokine P40, HP40, IL 9, IL9, IL-9, interleukin 9, P40, T cell growth factor P40
Accession Number	P15248
Source	Human Embryonic Kidney cells (HEK293). HEK293-derived IL-9 protein
Species Reactivity	human
Adventitious Virus	Master Cell Bank tested Negative for Adventitious Viruses

### Specifications


Test	Method	Specification
Activity	Dose-dependent stimulation of the proliferation of human MO7e cells (human megakaryoblastic leukemia cell line).	0.1-0.6 ng/mL
Molecular Mass	SDS-PAGE	38 to 45 kDa reduced, 32 to 40 kDa non-reduced, monomer, glycosylated
Purity	SDS-PAGE	>95%
Endotoxin	LAL	<0.1 EU/μg
Mycoplasma	PCR	Not Detected

### Activity Data



GMP Recombinant human IL-9 (HZ-1240-GMP) stimulates dose-dependent proliferation of the MO7e (human megakaryoblastic leukemia) cell line. Cell number was quantitatively assessed by PrestoBlue<sup>®</sup> Cell Viability Reagent. MO7e cells were treated with increasing concentrations of GMP recombinant IL-9 for 72 hours. The EC<sub>50</sub> was determined using a 4-parameter non-linear

### SDS-PAGE



Purity of recombinant human IL-9 was determined by SDS-polyacrylamide gel electrophoresis. The protein was resolved in an SDS-polyacrylamide gel in reducing and non-reducing conditions and stained using Coomassie blue.

Preparation	
Shipping Temperature	ambient temperature
Formulation	1x PBS, See Certificate of Analysis for details
Reconstitution	Briefly centrifuge the vial before opening. It is recommended to reconstitute the protein to 0.2 mg/mL in sterile 1x PBS pH 7.4 containing 0.1% endotoxin-free recombinant human serum albumin (HSA). Gently swirl or tap vial to mix.

Stability and Storage	Product Form	Temperature Conditions	Storage Time (From Date of Receipt)
	Lyophilized	-20°C to -80°C	Until Expiry Date
	Lyophilized	Room Temperature	2 weeks
	Reconstituted as per CofA	-20°C to -80°C	6 months
	Reconstituted as per CofA	4°C	1 week
Avoid repeated freeze-thaw cycles.			

## Proteintech GMP Quality Policy HumanKine® GMP Proteins

In vitro recombinant protein production can be prone to variability due to various reasons ranging from quality of raw materials to inconsistency in the process. Therefore, HumanKine®, a proteintech brand's GMP proteins are produced and tested under an ISO 13485 certified quality management system in a clean room facility. Proteintech manufactures the GMP HumanKine® products according to the applicable sections in the following documents: USP Chapter 1043 (Ancillary Materials for Cell, Gene, and Tissue-Engineered Products, USP Chapter 92 (Growth Factors and Cytokines Used in Cell Therapy Manufacturing), WHO TRS, No. 822, 1992 Annex 1 (Good Manufacturing Practices for Biological Products), Ph. Eur. General Chapter 5.2.12, and EudraLex – Volume 4 – Part IV (Guidelines on GMP specific to ATMPs). Proteintech strives to achieve the utmost quality GMP raw material ensuring all applicable guidelines are taken into consideration.

The QMS is built to provide our customers with consistent and pure product delivered by documented processes, qualified personnel, validated processes, qualified equipment, qualified vendors, and a stringent final product release process. Although the final product release process is important, Proteintech performs in-process QC steps after each major manufacturing stage. Production records and facilities may be available for an inspection by approved personnel.

Our GMP policy covers all the aspects of production; from raw materials, facility, equipment, and Instruments to training and personal hygiene of staff. It also guarantees that the process is explicit, validated and well documented for transparency and traceability.

[www.ptglab.com](http://www.ptglab.com)

Document #: FR-QA118-101 Rev 0  
Data Sheet Version #:

**Proteintech Group, Inc.**  
5500 Pearl Street, Suite 400 Rosemont, IL 60612  
t: 1-888-478-4522; f: 1-312-455-8408  
Email: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)