

IHC*easy* NAMPT Ready-To-Use IHC Kit

Catalog Number: **KHC0644**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Mouse Monoclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

Background

Nicotinamide phosphoribosyltransferase (NAMPT) has two usual synonyms termed Visfatin and PBEF. Its primary role is to catalyze the condensation of nicotinamide with 5-phosphoribosyl-1-pyrophosphate to yield nicotinamide mononucleotide, an intermediate in the biosynthesis of NAD, which is the rate limiting component in the mammalian NAD biosynthesis pathway. NAMPT is localized in cytoplasm and expressed in large amounts in bone marrow, liver tissue, and muscle tissues. NAMPT inhibits neutrophil apoptosis in experimental inflammation and clinical sepsis. NAMPT levels are altered in plasma of patients with type 2 diabetes mellitus (T2DM), and it is now evidenced that NAMPT may play a role in lipid metabolism.


Synonyms

1110035O14Rik, NAMPRtase, NAMPT, NAMPT/PBEF, PBEF, PBEF1, Pre B cell enhancing factor, VF, VISFATIN


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using KHC0644 (NAMPT/PBEF IHC Kit).




Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using KHC0644 (NAMPT/PBEF IHC Kit).




Immunohistochemical analysis of paraffin-embedded human appendicitis tissue slide using KHC0644 (NAMPT/PBEF IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using KHC0644 (NAMPT/PBEF IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using KHC0644 (NAMPT/PBEF IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC0644 (NAMPT/PBEF IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat lung tissue slide using KHC0644 (NAMPT/PBEF IHC Kit).