

IHC*easy* CORO1A Ready-To-Use IHC Kit

Catalog Number: **KHC0702**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	


Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Synonyms

CLABP, Clipin A, CLIPINA, CORO1, CORO1A, Coronin 1A, Coronin like protein A, Coronin like protein p57, HCORO1, p57, TACO


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human lymphoma tissue slide using KHC0702 (CORO1A IHC Kit).




Immunohistochemical analysis of paraffin-embedded human appendicitis tissue slide using KHC0702 (CORO1A IHC Kit).




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using KHC0702 (CORO1A IHC Kit).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using KHC0702 (CORO1A IHC Kit).



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC0702 (CORO1A IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse spleen tissue slide using KHC0702 (CORO1A IHC Kit).