

IHC*easy* MLKL Ready-To-Use IHC Kit

Catalog Number: **KHC2686**

General Information

Sample type:
FFPE tissue

Cited sample type:

Reactivity:
Human, Mouse, Rat

Cited Reactivity:

Assay type:
Immunohistochemistry

Primary antibody type:
Mouse Monoclonal

Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

Mixed lineage kinase domain like pseudokinase (MLKL), belongs to the protein kinase superfamily. MLKL plays a critical role in tumor necrosis factor (TNF)-induced necroptosis, a programmed cell death process, via interaction with receptor-interacting protein 3 (RIP3), which is a key signaling molecule in the necroptosis pathway. High levels of this protein and RIP3 are associated with inflammatory bowel disease in children.


Synonyms

MLKL,hMLKL,Mixed lineage kinase domain-like protei


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human rectal cancer tissue slide using KHC2686 (MLKL IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse spleen tissue slide using KHC2686 (MLKL IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse salivary gland tissue slide using KHC2686 (MLKL IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat lung tissue slide using KHC2686 (MLKL IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat salivary gland tissue slide using KHC2686 (MLKL IHC Kit).