

For Research Use Only

IFT140 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00029-3

Capture Antibody Information

Catalog Number:
82947-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230231A2

Reactivity:
Human

Immunogen Catalog Number:
Ag11119

Conjugate:
Unconjugated

Full name:
intraflagellar transport 140 homolog
(Chlamydomonas)

Gene ID:
9742

Detection Antibody Information

Catalog Number:
82947-6-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230231H2

Reactivity:
Human

GenBank:
BC035577

Immunogen Catalog Number:
Ag11119

Conjugate:
Unconjugated

Full name:
intraflagellar transport 140 homolog
(Chlamydomonas)

Gene ID:
9742

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00029-3 targets IFT140 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IFT140 Recombinant antibody, PBS Only (Capture) 82947-3-PBS (230231A2). 100 µg. Concentration 1 mg/ml.

Detection antibody: IFT140 Recombinant antibody, PBS Only (Detector) 82947-6-PBS (230231H2). 100 µg. Concentration 1 mg/ml.

Alternative IFT140 matched antibody pairs: MP00029-1, MP00029-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00029-3, IFT140 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82947-3-PBS. Detection antibody: 82947-6-PBS. Standard: Ag11119. Range: 2.5-20 ng/mL