

For Research Use Only

DDX60 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00038-2

Capture Antibody Information

Catalog Number:
83286-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240198H5

Reactivity:
Human

Immunogen Catalog Number:
Ag34626

Conjugate:
Unconjugated

Full name:
DEAD (Asp-Glu-Ala-Asp) box
polypeptide 60

Gene ID:
55601

Detection Antibody Information

Catalog Number:
83286-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240198B3

Reactivity:
Human

GenBank:
BC038115

Immunogen Catalog Number:
Ag34626

Conjugate:
Unconjugated

Full name:
DEAD (Asp-Glu-Ala-Asp) box
polypeptide 60

Gene ID:
55601

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00038-2 targets DDX60 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDX60 Recombinant antibody, PBS Only (Capture) 83286-2-PBS (240198H5). 100 µg. Concentration 1 mg/ml.

Detection antibody: DDX60 Recombinant antibody, PBS Only (Detector) 83286-1-PBS (240198B3). 100 µg. Concentration 1 mg/ml.

Alternative DDX60 matched antibody pairs: MP00038-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00038-2, DDX60 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83286-2-PBS. Detection antibody: 83286-1-PBS. Standard: Ag34626. Range: 0.78-100ng/mL.