

For Research Use Only

CAMP Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00042-2

Capture Antibody Information

Catalog Number:
82940-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230157E9

Reactivity:
Human

Immunogen Catalog Number:
Ag2622

Conjugate:
Unconjugated

Full name:
cathelicidin antimicrobial peptide

Gene ID:
820

Detection Antibody Information

Catalog Number:
82940-6-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230157H2

Reactivity:
Human

GenBank:
BC055089

Immunogen Catalog Number:
Ag2622

Conjugate:
Unconjugated

Full name:
cathelicidin antimicrobial peptide

Gene ID:
820

Applications

Tested Applications:
Sandwich ELISA

Range:
0.39-25 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00042-2 targets CAMP in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: CAMP Recombinant antibody, PBS Only (Capture) 82940-4-PBS (230157E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: CAMP Recombinant antibody, PBS Only (Detector) 82940-6-PBS (230157H2). 100 µg. Concentration 1 mg/mL.

Alternative CAMP matched antibody pairs: MP00042-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00042-2, CAMP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82940-4-PBS. Detection antibody: 82940-6-PBS. Standard: Ag2622. Range: 0.39-25 ng/mL.