

For Research Use Only

# STT3A Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00080-1

## Capture Antibody Information

**Catalog Number:**  
83077-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230426G7  
**Reactivity:**  
Human  
**Immunogen Catalog Number:**  
Ag31994

**Conjugate:**  
Unconjugated  
**Full name:**  
STT3, subunit of the oligosaccharyltransferase complex, homolog A (*S. cerevisiae*)  
**Gene ID:**  
3703

## Detection Antibody Information

**Catalog Number:**  
83077-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230426D2  
**Reactivity:**  
Human  
**GenBank:**  
BC020965  
**Immunogen Catalog Number:**  
Ag31994

**Conjugate:**  
Unconjugated  
**Full name:**  
STT3, subunit of the oligosaccharyltransferase complex, homolog A (*S. cerevisiae*)  
**Gene ID:**  
3703

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00080-1 targets STT3A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STT3A Recombinant antibody, PBS Only (Capture) 83077-3-PBS (230426G7). 100 µg. Concentration 1 mg/mL.

Detection antibody: STT3A Recombinant antibody, PBS Only (Detector) 83077-2-PBS (230426D2). 100 µg. Concentration 1 mg/mL.

Alternative STT3A matched antibody pairs: MP00080-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00080-1, STT3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83077-3-PBS. Detection antibody: 83077-2-PBS. Standard: Ag31994. Range: 0.625-20 ng/mL