

For Research Use Only

# UBN1-Specific Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00095-3

## Capture Antibody Information

**Catalog Number:**  
83001-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240014H3  
**Reactivity:**  
Human

**Conjugate:**  
Unconjugated  
**Full name:**  
ubিনuclein 1  
**Gene ID:**  
29855

## Detection Antibody Information

**Catalog Number:**  
83001-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240014A8  
**Reactivity:**  
Human  
**GenBank:**  
NM\_016936

**Conjugate:**  
Unconjugated  
**Full name:**  
ubিনuclein 1  
**Gene ID:**  
29855

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.078-5 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00095-3 targets UBN1-Specific in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: UBN1-Specific Recombinant antibody, PBS Only (Capture) 83001-3-PBS (240014H3). 100 µg. Concentration 1 mg/ml.

Detection antibody: UBN1-Specific Recombinant antibody, PBS Only (Detector) 83001-1-PBS (240014A8). 100 µg. Concentration 1 mg/ml.

Alternative UBN1-Specific matched antibody pairs: MP00095-1, MP00095-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00095-3, UBN1-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83001-3-PBS. Detection antibody: 83001-1-PBS. Standard: SY00809. Range: 0.078-5 ng/mL.