

For Research Use Only

# DDX39 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00111-1

## Capture Antibody Information

**Catalog Number:**  
83083-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230375G6  
**Reactivity:**  
Human  
**Immunogen Catalog Number:**  
Ag2311

**Conjugate:**  
Unconjugated  
**Full name:**  
DEAD (Asp-Glu-Ala-Asp) box polypeptide 39  
**Gene ID:**  
10212

## Detection Antibody Information

**Catalog Number:**  
83083-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230375E8  
**Reactivity:**  
Human  
**GenBank:**  
BC032128  
**Immunogen Catalog Number:**  
Ag2311

**Conjugate:**  
Unconjugated  
**Full name:**  
DEAD (Asp-Glu-Ala-Asp) box polypeptide 39  
**Gene ID:**  
10212

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00111-1 targets DDX39 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDX39 Recombinant antibody, PBS Only (Capture) 83083-3-PBS (230375G6). 100 µg. Concentration 1 mg/ml.

Detection antibody: DDX39 Recombinant antibody, PBS Only (Detector) 83083-2-PBS (230375E8). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00111-1, DDX39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83083-3-PBS. Detection antibody: 83083-2-PBS. Standard: Ag2311. Range: 0.625-20 ng/mL