

For Research Use Only

Mouse KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00116-1

Capture Antibody Information

Catalog Number:
83064-6-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230336H6

Reactivity:
mouse

Conjugate:
Unconjugated

Full name:
hepatitis A virus cellular receptor 1

Gene ID:
171283

Detection Antibody Information

Catalog Number:
83064-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230336D2

Reactivity:
mouse

GenBank:
NM_134248

Conjugate:
Unconjugated

Full name:
hepatitis A virus cellular receptor 1

Gene ID:
171283

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00116-1 targets KIM-1/HAVCR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse KIM-1/HAVCR1 Recombinant antibody, PBS Only (Capture) 83064-6-PBS (230336H6). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse KIM-1/HAVCR1 Recombinant antibody, PBS Only (Detector) 83064-3-PBS (230336D2). 100 µg. Concentration 1 mg/ml.

Alternative KIM-1/HAVCR1 matched antibody pairs: MP00116-2, MP00116-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00116-1, Mouse KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83064-6-PBS. Detection antibody: 83064-3-PBS. Standard: Eg0571. Range: 0.625-10 ng/mL.