

For Research Use Only



NFX1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00127-3

Capture Antibody Information

Catalog Number:

82924-5-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230211G5

Reactivity:

Human

Immunogen Catalog Number:

Ag9448

Conjugate:

Unconjugated

Full name:

nuclear transcription factor, X-box binding 1

Gene ID:

4799

Detection Antibody Information

Catalog Number:

82924-6-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230211H6

Reactivity:

Human

GenBank:

BC012151

Immunogen Catalog Number:

Ag9448

Conjugate:

Unconjugated

Full name:

nuclear transcription factor, X-box binding 1

Gene ID:

4799

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00127-3 targets NFX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFX1 Recombinant antibody, PBS Only (Capture) 82924-5-PBS (230211G5). 100 µg. Concentration 1 mg/ml.

Detection antibody: NFX1 Recombinant antibody, PBS Only (Detector) 82924-6-PBS (230211H6). 100 µg. Concentration 1 mg/ml.

Alternative NFX1 matched antibody pairs: MP00127-1, MP00127-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00127-3, NFX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82924-5-PBS. Detection antibody: 82924-6-PBS. Standard: Ag9448. Range: 0.625-20 ng/mL.