

For Research Use Only

DDX19A Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00140-1

Capture Antibody Information

Catalog Number:
83106-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230232D3
Reactivity:
Human
Immunogen Catalog Number:
Ag7362

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-As) box polypeptide 19A
Gene ID:
55308

Detection Antibody Information

Catalog Number:
83106-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230232F11
Reactivity:
Human
GenBank:
BC005162
Immunogen Catalog Number:
Ag7362

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-As) box polypeptide 19A
Gene ID:
55308

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00140-1 targets DDX19A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDX19A Recombinant antibody, PBS Only (Capture) 83106-3-PBS (230232D3). 100 µg. Concentration 1 mg/ml.

Detection antibody: DDX19A Recombinant antibody, PBS Only (Detector) 83106-6-PBS (230232F11). 100 µg. Concentration 1 mg/ml.

Alternative DDX19A matched antibody pairs: MP00140-2, MP00140-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00140-1, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-3-PBS. Detection antibody: 83106-6-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL.