

For Research Use Only

CBX5 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00181-3

Capture Antibody Information

Catalog Number:

83258-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240145G5

Reactivity:

Human

Immunogen Catalog Number:

Ag34810

Conjugate:

Unconjugated

Full name:

chromobox homolog 5 (HP1 alpha homolog, Drosophila)

Gene ID:

23468

Detection Antibody Information

Catalog Number:

83258-5-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240145F11

Reactivity:

Human

GenBank:

BC006821

Immunogen Catalog Number:

Ag34810

Conjugate:

Unconjugated

Full name:

chromobox homolog 5 (HP1 alpha homolog, Drosophila)

Gene ID:

23468

Applications

Tested Applications:

Cytometric bead array

Range:

3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00181-3 targets CBX5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CBX5 Recombinant antibody, PBS Only (Capture) 83258-2-PBS (240145G5). 100 µg. Concentration 1 mg/ml.

Detection antibody: CBX5 Recombinant antibody, PBS Only (Detector) 83258-5-PBS (240145F11). 100 µg. Concentration 1 mg/ml.

Alternative CBX5 matched antibody pairs: MP00181-1, MP00181-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00181-3, CBX5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83258-2-PBS. Detection antibody: 83258-5-PBS. Standard: Ag34810. Range: 3.125-100 ng/mL.