

For Research Use Only

NRG1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00232-6

Capture Antibody Information

Catalog Number:
83251-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240067D3

Reactivity:
human

Immunogen Catalog Number:
Ag35114

Conjugate:
Unconjugated

Full name:
neuregulin 1

Gene ID:
3084

Detection Antibody Information

Catalog Number:
83251-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240067A4

Reactivity:
human

GenBank:
BC007675

Immunogen Catalog Number:
Ag35114

Conjugate:
Unconjugated

Full name:
neuregulin 1

Gene ID:
3084

Applications

Tested Applications:
Sandwich ELISA

Range:
78.1-5000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00232-6 targets NRG1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: NRG1 Recombinant antibody, PBS Only (Capture) 83251-4-PBS (240067D3). 100 µg. Concentration 1 mg/mL.

Detection antibody: NRG1, isoform Alpha Recombinant antibody, PBS Only (Detector) 83251-1-PBS (240067A4). 100 µg. Concentration 1 mg/mL.

Alternative NRG1 matched antibody pairs: MP00232-1, MP00232-2, MP00232-3, MP00232-4, MP00232-5

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00232-6, Human NRG1 Recombinant Matched Antibody Pair - PBS only. 83251-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag35114. 83251-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean NRG1 concentration was determined to be 2,318.1 pg/mL in A549 cell extract based on a 2.5 mg/mL extract load.