

For Research Use Only



PPARG Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00256-1

Capture Antibody Information

Catalog Number:
81490-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230374C2
Reactivity:
Human
Immunogen Catalog Number:
Ag10005

Conjugate:
Unconjugated
Full name:
peroxisome proliferator-activated receptor gamma
Gene ID:
5468

Detection Antibody Information

Catalog Number:
81490-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230374E7
Reactivity:
Human
GenBank:
BC006811
Immunogen Catalog Number:
Ag10005

Conjugate:
Unconjugated
Full name:
peroxisome proliferator-activated receptor gamma
Gene ID:
5468

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00256-1 targets PPARG in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PPARG Recombinant antibody, PBS Only (Capture) 81490-3-PBS (230374C2). 100 µg. Concentration 1 mg/ml.

Detection antibody: PPARG Recombinant antibody, PBS Only (Detector) 81490-4-PBS (230374E7). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00256-1, PPAR γ Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81490-3-PBS. Detection antibody: 81490-4-PBS. Standard: Ag10005. Range: 1.25-20 ng/mL