

For Research Use Only

# TPR Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00258-2

## Capture Antibody Information

**Catalog Number:**  
83177-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230413C9  
**Reactivity:**  
Human

**Conjugate:**  
Unconjugated  
**Full name:**  
translocated promoter region (to activated MET oncogene)  
**Gene ID:**  
7175

## Detection Antibody Information

**Catalog Number:**  
83177-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230413C3  
**Reactivity:**  
Human  
**GenBank:**  
NM\_003292

**Conjugate:**  
Unconjugated  
**Full name:**  
translocated promoter region (to activated MET oncogene)  
**Gene ID:**  
7175

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.312-5 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00258-2 targets TPR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TPR Recombinant antibody, PBS Only (Capture) 83177-4-PBS (230413C9). 100 µg. Concentration 1 mg/mL.

Detection antibody: TPR Recombinant antibody, PBS Only (Detector) 83177-3-PBS (230413C3). 100 µg. Concentration 1 mg/mL.

Alternative TPR matched antibody pairs: MP00258-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00258-2, TPR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83177-4-PBS. Detection antibody: 83177-3-PBS. Standard: SY01141. Range: 0.312-5 ng/mL