

For Research Use Only

# TGS1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00299-1

## Capture Antibody Information

Catalog Number:  
83140-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
230427F6

Reactivity:  
Human

Immunogen Catalog Number:  
Ag3683

Conjugate:  
Unconjugated

Full name:  
trimethylguanosine synthase  
homolog (*S. cerevisiae*)

Gene ID:  
96764

## Detection Antibody Information

Catalog Number:  
83140-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
230427G7

Reactivity:  
Human

GenBank:  
BC011999

Immunogen Catalog Number:  
Ag3683

Conjugate:  
Unconjugated

Full name:  
trimethylguanosine synthase  
homolog (*S. cerevisiae*)

Gene ID:  
96764

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00299-1 targets TGS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TGS1 Recombinant antibody, PBS Only (Capture) 83140-2-PBS (230427F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: TGS1 Recombinant antibody, PBS Only (Detector) 83140-4-PBS (230427G7). 100 µg. Concentration 1 mg/mL.

Alternative TGS1 matched antibody pairs: MP00299-2, MP00299-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00299-1, TGS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83140-2-PBS. Detection antibody: 83140-4-PBS. Standard: Ag3683. Range: 0.625-20 ng/mL