

For Research Use Only

Mouse TREM-1/CD354 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00316-1

Capture Antibody Information

Catalog Number:
83294-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240219D3
Reactivity:
mouse

Conjugate:
Unconjugated
Full name:
triggering receptor expressed on myeloid cells 1
Gene ID:
58217

Detection Antibody Information

Catalog Number:
83294-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240219H1
Reactivity:
mouse
GenBank:
NM_021406.5

Conjugate:
Unconjugated
Full name:
triggering receptor expressed on myeloid cells 1
Gene ID:
58217

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00316-1 targets TREM-1/CD354 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse TREM-1/CD354 Recombinant antibody, PBS Only (Capture) 83294-1-PBS (240219D3). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse TREM-1/CD354 Recombinant antibody, PBS Only (Detector) 83294-2-PBS (240219H1). 100 µg. Concentration 1 mg/ml.

Alternative TREM-1/CD354 matched antibody pairs: MP00988-1, MP00988-2, MP00988-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00316-1, MOUSE Trem1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83294-1-PBS. Detection antibody: 83294-2-PBS. Standard: Eg0874. Range: 1.25-40 ng/mL