

For Research Use Only

# IGFBP7 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00319-2

## Capture Antibody Information

**Catalog Number:**  
83292-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240228F1  
**Reactivity:**  
Human

**Conjugate:**  
Unconjugated  
**Full name:**  
insulin-like growth factor binding protein 7  
**Gene ID:**  
3490

## Detection Antibody Information

**Catalog Number:**  
83292-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240228C4  
**Reactivity:**  
Human  
**GenBank:**  
NM\_001553.2

**Conjugate:**  
Unconjugated  
**Full name:**  
insulin-like growth factor binding protein 7  
**Gene ID:**  
3490

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.25-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00319-2 targets IGFBP7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IGFBP7 Recombinant antibody, PBS Only (Capture) 83292-3-PBS (240228F1). 100 µg. Concentration 1 mg/ml.

Detection antibody: IGFBP7 Recombinant antibody, PBS Only (Detector) 83292-1-PBS (240228C4). 100 µg. Concentration 1 mg/ml.

Alternative IGFBP7 matched antibody pairs: MP00319-1, MP00319-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00319-2, IGFBP7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83292-3-PBS. Detection antibody: 83292-1-PBS. Standard: Eg0815. Range: 1.25-40 ng/mL