For Research Use Only

NUPR1 Recombinant Matched Antibody Pair, PBS Only

Proteintech®
Antibodies | ELISA kits | Proteins
www.ptglab.com

Catalog Number: MP00339-1

Capture Antibody Information Catalog Number: Clone ID: 83315-1-PBS 240149A8
Host: Reactivity:

Isotype: Immunogen Catalog Number: Gene ID: IgG Ag7039 Gene ID: 26471

Human

Purification Method: Protein A purification

Rabbit

Detection Antibody Information

Catalog Number: Clone ID: Conjugate: 83315-3-PBS 240149F4 Unconjugated Host: Reactivity: Full name: Rabbit Human nuclear protein 1 Isotype: GenBank: Gene ID: IgG BC002434 26471

Purification Method: Immunogen Catalog Number:

Protein A purification Ag7039

Applications

Tested Applications: Rang

Cytometric bead array 3.13-100 ng/mL (Cytometric Bead

Array)

Recommended Dilutions:

Conjugate:

Full name:

Unconjugated

nuclear protein 1

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00339-1 targets NUPR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NUPR1 Recombinant antibody, PBS Only (Capture) 83315-1-PBS (240149A8). 100 µg. Concentration 1 mgl/ml.

Detection antibody: NUPR1 Recombinant antibody, PBS Only (Detector) 83315-3-PBS (240149F4). 100 µg. Concentration 1 mgl/ml.

Alternative NUPR1 matched antibody pairs: MP00339-2, MP00339-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage: Store at -80°C. Storage buffer: PBS only

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com W: ptglab.com

Selected Validation Data



Cytometric bead array standard curve of MP00339-1, NUPR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83315-1-PBS. Detection antibody: 83315-3-PBS. Standard: Ag7039. Range: 3.13-100 ng/mL