

For Research Use Only

NLRC5 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00402-3

Capture Antibody Information

Catalog Number:
83384-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240338D8

Reactivity:
Human

Immunogen Catalog Number:
Ag34508

Conjugate:
Unconjugated

Full name:
NLR family, CARD domain containing 5

Gene ID:
84166

Detection Antibody Information

Catalog Number:
83384-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240338B12

Reactivity:
Human

GenBank:
BC063566

Immunogen Catalog Number:
Ag34508

Conjugate:
Unconjugated

Full name:
NLR family, CARD domain containing 5

Gene ID:
84166

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00402-3 targets NLRC5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NLRC5 Recombinant antibody, PBS Only (Capture) 83384-2-PBS (240338D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: NLRC5 Recombinant antibody, PBS Only (Capture/Detector) 83384-3-PBS (240338B12). 100 µg. Concentration 1 mg/ml.

Alternative NLRC5 matched antibody pairs: MP00402-1, MP00402-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00402-3, NLRC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83384-2-PBS. Detection antibody: 83384-3-PBS. Standard: Ag34508. Range: 1.56-100 ng/mL