

For Research Use Only

GNB3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00443-2

Capture Antibody Information

Catalog Number:
83452-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240281C5

Reactivity:
Human

Immunogen Catalog Number:
Ag0127

Conjugate:
Unconjugated

Full name:
guanine nucleotide binding protein (G protein), beta polypeptide 3

Gene ID:
2784

Detection Antibody Information

Catalog Number:
83452-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240281B5

Reactivity:
Human

GenBank:
BC000115

Immunogen Catalog Number:
Ag0127

Conjugate:
Unconjugated

Full name:
guanine nucleotide binding protein (G protein), beta polypeptide 3

Gene ID:
2784

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00443-2 targets GNB3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GNB3 Recombinant antibody, PBS Only (Capture) 83452-3-PBS (240281C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: GNB3 Recombinant antibody, PBS Only (Capture/Detector) 83452-1-PBS (240281B5). 100 µg. Concentration 1 mg/mL.

Alternative GNB3 matched antibody pairs: MP00443-1, MP00443-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00443-2, GNB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83452-3-PBS. Detection antibody: 83452-1-PBS. Standard: Ag0127. Range: 0.78-100 ng/mL