

For Research Use Only



IFITM2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00456-3

Capture Antibody Information

Catalog Number:
83455-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240318F6

Reactivity:
Human

Immunogen Catalog Number:
Ag3451

Conjugate:
Unconjugated

Full name:
interferon induced transmembrane protein 2 (1-8D)

Gene ID:
10581

Detection Antibody Information

Catalog Number:
83455-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240318B7

Reactivity:
Human

GenBank:
BC009696

Immunogen Catalog Number:
Ag3451

Conjugate:
Unconjugated

Full name:
interferon induced transmembrane protein 2 (1-8D)

Gene ID:
10581

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00456-3 targets IFITM2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IFITM2 Recombinant antibody, PBS Only (Capture/Detector) 83455-4-PBS (240318F6). 100 µg. Concentration 1 mg/ml.

Detection antibody: IFITM2 Recombinant antibody, PBS Only (Detector) 83455-1-PBS (240318B7). 100 µg. Concentration 1 mg/ml.

Alternative IFITM2 matched antibody pairs: MP00456-1, MP00456-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00456-3, IFITM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83455-4-PBS. Detection antibody: 83455-1-PBS. Standard: Ag3451. Range: 0.78-100 ng/mL