

For Research Use Only

# PMS1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00457-1

## Capture Antibody Information

**Catalog Number:**  
83456-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240408D8  
**Reactivity:**  
Human  
**Immunogen Catalog Number:**  
Ag28132

**Conjugate:**  
Unconjugated  
**Full name:**  
PMS1 postmeiotic segregation increased 1 (*S. cerevisiae*)  
**Gene ID:**  
5378

## Detection Antibody Information

**Catalog Number:**  
83456-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240408H3  
**Reactivity:**  
Human  
**GenBank:**  
BC096331  
**Immunogen Catalog Number:**  
Ag28132

**Conjugate:**  
Unconjugated  
**Full name:**  
PMS1 postmeiotic segregation increased 1 (*S. cerevisiae*)  
**Gene ID:**  
5378

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.156-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00457-1 targets PMS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PMS1 Recombinant antibody, PBS Only (Capture) 83456-3-PBS (240408D8). 100 µg. Concentration 1 mg/mL.

Detection antibody: PMS1 Recombinant antibody, PBS Only (Detector) 83456-4-PBS (240408H3). 100 µg. Concentration 1 mg/mL.

Alternative PMS1 matched antibody pairs: MP00457-2, MP00457-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00457-1, PMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83456-3-PBS. Detection antibody: 83456-4-PBS. Standard: Ag28132. Range: 0.156-20 ng/mL