

For Research Use Only



TPSAB1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00516-1

Capture Antibody Information

Catalog Number:
83518-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240432E12

Reactivity:
Human

Immunogen Catalog Number:
Ag4167

Conjugate:
Unconjugated

Full name:
tryptase alpha/beta 1

Gene ID:
7177

Detection Antibody Information

Catalog Number:
83518-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240432A2

Reactivity:
Human

GenBank:
BC028059

Immunogen Catalog Number:
Ag4167

Conjugate:
Unconjugated

Full name:
tryptase alpha/beta 1

Gene ID:
7177

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00516-1 targets TPSAB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TPSAB1 Recombinant antibody, PBS Only (Capture) 83518-2-PBS (240432E12). 100 µg. Concentration 1 mg/ml.

Detection antibody: TPSAB1 Recombinant antibody, PBS Only (Detector) 83518-1-PBS (240432A2). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00516-1, TPSAB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83518-2-PBS. Detection antibody: 83518-1-PBS. Standard: Ag4167. Range: 0.78-100 ng/mL.