

For Research Use Only

Mouse ICAM-1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00543-3

Capture Antibody Information

Catalog Number:

82827-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240249A10

Reactivity:

Mouse

Conjugate:

Unconjugated

Full name:

intercellular adhesion molecule 1

Gene ID:

15894

Detection Antibody Information

Catalog Number:

82827-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240249G1

Reactivity:

Mouse

GenBank:

NM_010493

Conjugate:

Unconjugated

Full name:

intercellular adhesion molecule 1

Gene ID:

15894

Applications

Tested Applications:

Cytometric bead array

Range:

0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00543-3 targets ICAM-1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse ICAM1 Recombinant antibody, PBS Only (Capture) 82827-2-PBS (240249A10). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse ICAM1 Recombinant antibody, PBS Only (Detector) 82827-4-PBS (240249G1). 100 µg. Concentration 1 mg/ml.

Alternative ICAM-1 matched antibody pairs: MP00028-1, MP00028-2, MP00028-3, MP00028-4, MP00028-5, MP00543-1, MP00543-2, MP50085-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00543-3, MOUSE ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82827-2-PBS. Detection antibody: 82827-4-PBS. Standard: Eg0482. Range: 0.313-40 ng/mL.